


Anti-Ldb3 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - O75112

Item Number: CM0027868



Introduction

High-quality Anti-LDB3 Rabbit Polyclonal Antibody (O75112) validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat LDB3. Ideal for striated muscle, cardiomyopathy, and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	LDB3 Rabbit pAb
Remarks/Alias	MF4, ZASP, CMD1C, CMH24, CMPD3, LVNC3, CYPHER, LDB3Z1, LDB3Z4, ORACLE, PDLIM6, LDB3
Species	Human
GeneID (Human)	11155
GeneID	11155
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human LDB3 (NP_001073585.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75112
Protein Weight	77kDa
Shipping	Ice bag