


Anti-Ldha Monoclonal Antibody For Wb And Elisa - P00338

Item Number: CM0012008



Introduction

Rabbit monoclonal antibody targeting human LDHA, validated for Western blotting and ELISA. Detects LDHA in human, mouse, and rat samples; ideal for glycolysis and metabolic disorder research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] LDHA Rabbit mAb
Remarks/Alias	LDHM; GSD11; PIG19; HEL-S-133P
Species	Human
Gene ID (Human)	3939
Gene ID	3939
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LDHA (NP_005557.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00338
Protein Weight	26-39 kDa
Shipping	Ice bag