


Anti-Ldha Monoclonal Antibody For Wb And Elisa - P00338

Item Number: CM0007143



Introduction

Rabbit monoclonal antibody targeting human lactate dehydrogenase A (LDHA), validated in WB and ELISA. Detects human, mouse, and rat LDHA. Ideal for metabolism, glycolysis, and cancer biomarker research.

[Learn More](#)

Parameter	Value
Product Name	LDHA Rabbit mAb
Remarks/Alias	LDHM; GSD11; PIG19; HEL-S-133P; LDHA
Species	Human
GeneID (Human)	3939
GeneID	3939
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00338
Protein Weight	26-39 kDa
Shipping	Ice bag