


Anti-Ldha Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P00338

Item Number: CM0014338



Introduction

Rabbit polyclonal antibody against human LDHA. Validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Ideal for studying glycolysis, hypoxia, and cancer metabolism.

[Learn More](#)

Parameter	Value
Product Name	LDHA Rabbit pAb
Remarks/Alias	LDHM; GSD11; PIG19; HEL-S-133P; LDHA
Species	Human
Gene ID (Human)	3939
Gene ID	3939
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-332 of human LDHA (NP_005557.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00338
Protein Weight	26-39 kDa
Shipping	Ice bag