


Anti-Ldha Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P00338

Item Number: CM0011551



Introduction

Rabbit monoclonal antibody against human LDHA, validated by knockdown for Western blot, IF/ICC, and ELISA. Detects endogenous LDHA in human, mouse, and rat samples. Predicted molecular weight 26-39 kDa. Ideal for cancer metabolism and glycolysis research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] LDHA Rabbit mAb
Remarks/Alias	LDHM; GSD11; PIG19; HEL-S-133P
Species	Human
Gene ID (Human)	3939
Gene ID	3939
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-332 of human LDHA (NP_005557.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00338
Protein Weight	26-39 kDa
Shipping	Ice bag