


Anti-Ldha Rabbit Polyclonal Antibody For Ihc-P/If/Icc/Elisa - P00338

Item Number: CM0019176



Introduction

Rabbit polyclonal anti-LDHA antibody validated for IHC-P, IF/ICC, and ELISA. Detects human and mouse LDHA, a glycolysis enzyme linked to cancer metabolism and glycogen storage disease.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] LDHA Rabbit pAb
Remarks/Alias	LDHM; GSD11; PIG19; HEL-S-133P; LDHA
Species	Human
Gene ID (Human)	3939
Gene ID	3939
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human LDHA (NP_001158886.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00338
Protein Weight	26-39 kDa
Shipping	Ice bag