

Anti-Ldhal6B Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9Byz2

Item Number: CM0018460

Introduction

Rabbit polyclonal antibody targeting human LDHAL6B, validated for WB, IF-P, ELISA. Testis-specific lactate dehydrogenase catalyzing lactate-pyruvate interconversion. Cross-reacts with mouse and rat. Ideal for metabolism and fertility research.

[Learn More](#)

Parameter	Value
Product Name	LDHAL6B Rabbit pAb
Remarks/Alias	LDHL; LDH6B; LDHAL6; LDHAL6B
Species	Human
GeneID (Human)	92483
GeneID	92483
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human LDHAL6B (NP_149972.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9BYZ2
Protein Weight	42kDa
Shipping	Ice bag