


Anti-Ldhb Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P07195

Item Number: CM0025333



Introduction

Rabbit polyclonal antibody targeting human LDHB, suitable for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. LDHB is a key enzyme in anaerobic glycolysis and a cardiac tissue marker.

[Learn More](#)

Parameter	Value
Product Name	LDHB Rabbit pAb
Remarks/Alias	LDH-B; LDH-H; LDHBD; TRG-5; HEL-S-281; LDHB
Species	Human
GeneID (Human)	3945
GeneID	3945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human LDHB (NP_002291.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07195
Protein Weight	37kDa
Shipping	Ice bag