


Anti-Ldhc Polyclonal Antibody For Wb, If-P, Elisa - P07864

Item Number: CM0014161



Introduction

Anti-LDHC Polyclonal Antibody for WB, IF-P, ELISA applications; cross-reacts with Human, Mouse, Rat. Targets L-lactate dehydrogenase C chain, a sperm motility-related protein and cancer/testis antigen 32.

[Learn More](#)

Parameter	Value
Product Name	LDHC Rabbit pAb
Remarks/Alias	CT32; LDH3; LDHX; LDHC
Species	Human
Gene ID	3948
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 150-250 of human LDHC (NP_002292.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07864
Protein Weight	36kDa
Shipping	Ice bag