

Anti-Ldl Receptor (Ldlr) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01130

Item Number: CM0016883



Introduction

Rabbit polyclonal antibody targeting human LDL Receptor (LDLR), validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat LDLR. Ideal for cholesterol metabolism and viral infection research.

[Learn More](#)

Parameter	Value
Product Name	LDL Receptor (LDLR) Rabbit pAb
Remarks/Alias	FH; FHC; FHCL1; LDLCQ2; LDL Receptor (LDLR)
Species	Human
GeneID (Human)	3949
GeneID	3949
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 761-860 of human LDL Receptor (LDLR) (NP_000518.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01130
Protein Weight	95kDa
Shipping	Ice bag