

Anti-Ldl Receptor (Ldlr) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01130

Item Number: CM0020740



Introduction

High-quality rabbit polyclonal antibody against human LDL Receptor (LDLR), validated for WB, IHC-P, IF/ICC, and ELISA, detecting human, mouse, and rat samples. Ideal for cholesterol metabolism, viral entry studies, and cardiovascular disease research.

[Learn More](#)

Parameter	Value
Product Name	LDL Receptor (LDLR) Rabbit pAb
Remarks/Alias	FH; FHC; FHCL1; LDLCQ2; LDL Receptor (LDLR)
Species	Human
GeneID	3949 (Human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 761-860 of human LDL Receptor (LDLR) (NP_000518.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01130
Protein Weight	95kDa
Shipping	Ice bag