


Anti-Ldlrap1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Sw96

Item Number: CM0026144



Introduction

Rabbit polyclonal antibody against human LDLRAP1 (ARH), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studies of LDLR endocytosis, hypercholesterolemia, and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	LDLRAP1 Rabbit pAb
Remarks/Alias	ARH; ARH1; ARH2; FHCB1; FHCB2; FHCL4; LDLRAP1
Species	Human
GeneID	26119
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-308 of human LDLRAP1 (NP_056442.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5SW96
Protein Weight	34kDa
Shipping	Ice bag