

Anti-Lef1 Polyclonal Antibody For Wb, Ip, If-P, Elisa - Q9Uju2

Item Number: CM0018202

Introduction

LEF1 Rabbit pAb (polyclonal) validated in WB, IP, IF-P, ELISA. Detects LEF1 (44 kDa) in human, mouse, rat. Recombinant immunogen (aa 20-132). Ideal for Wnt/transcription research.

[Learn More](#)

Parameter	Value
Product Name	LEF1 Rabbit pAb
Remarks/Alias	LEF-1; TCF10; TCF7L3; TCF1ALPHA; LEF1
Species	Human
GeneID (human)	51176
GeneID	51176
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-132 of human LEF1 (NP_057353.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJU2
Protein Weight	44 kDa
Shipping	Ice bag