

Anti-Lef1 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Uju2

Item Number: CM0007843

Introduction

Rabbit monoclonal LEF1 antibody validated for WB and ELISA, cross-reactive with human and mouse. Ideal for studying Wnt signaling and cancer biology in thymus, colon, and pancreatic tissues.

[Learn More](#)

Parameter	Value
Product Name	LEF1 Rabbit mAb
Remarks/Alias	LEF-1; TCF10; TCF7L3; TCF1ALPHA; LEF1
Species	Human
GeneID (Human)	51176
GeneID	51176
Immunogen	Recombinant protein of a sequence within amino acids 50-150 of human LEF1 (Q9UJU2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UJU2
Protein Weight	44kDa
Shipping	Ice bag