

Anti-Lef1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9Uju2

Item Number: CM0007591

Introduction

High-performance anti-LEF1 rabbit monoclonal antibody validated for WB, IP, and ELISA in human and mouse samples. Suitable for Wnt signaling, cancer, and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	LEF1 Rabbit mAb
Remarks/Alias	LEF-1; TCF10; TCF7L3; TCF1ALPHA; LEF1
Species	Human
Gene ID (Human)	51176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human LEF1 (NP_057353.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UJU2
Protein Weight	44kDa
Shipping	Ice bag