

Anti-Lefty1 Rabbit Polyclonal Antibody For Wb And Elisa - O75610

Item Number: CM0021079

Introduction

High-quality rabbit polyclonal antibody against human LEFTY1, validated for Western blot and ELISA. Detects human and mouse left-right determination factor 1, a key developmental regulator.

[Learn More](#)

Parameter	Value
Product Name	LEFTY1 Rabbit pAb
Remarks/Alias	LEFTB; LEFTYB; LEFTY1
Species	Human
GeneID (Human)	10637
GeneID	10637
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-366 of human LEFTY1 (NP_066277.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75610
Protein Weight	41kDa
Shipping	Ice bag