

Anti-Leghemoglobin (Mtlbs) Rabbit Polyclonal Antibody For Wb And Elisa - P27992/P27993

Item Number: CM0019443

Introduction

Rabbit polyclonal antibody targeting Medicago truncatula leghemoglobin (MtLBs/LB2), suitable for WB and ELISA. Recognizes MtLb2 and MtLb6 (P27992/P27993). Ideal for nitrogen fixation and nodulation studies.

[Learn More](#)

Parameter	Specification
Product Name	MtLBs Rabbit pAb
Remarks/Alias	LB2; MTlb2; MtLBs
Species	Medicago truncatula
Gene ID	11433347
Immunogen	Synthetic Peptide Two synthetic peptides corresponding to the sequence within amino acids 7-21 and 100-113 of Medicago truncatula Leghemoglobin 2 (XP_003588689.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Medicago truncatula
SWISS	P27992/P27993
Protein Weight	15kDa
Shipping	Ice bag