


Anti-Legumain (Lgmn) Monoclonal Antibody For Wb, Elisa - Q99538

Item Number: CM0008180



Introduction

Rabbit monoclonal antibody against Legumain (LGMN), validated for Western blot and ELISA. Detects human, mouse, rat; recombinant immunogen (aa 210-302). Suitable for studying protease function in lysosomal protein degradation and antigen presentation.

[Learn More](#)

Parameter	Value
Product Name	Legumain (LGMN) Rabbit mAb
Remarks/Alias	AEP; LGMN1; PRSC1; Legumain (LGMN)
Species	Human
GeneID (Human)	5641
GeneID	5641
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-302 of human Legumain (LGMN) (Q99538)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99538
Protein Weight	49 kDa
Shipping	Ice bag