


Anti-Legumain (Lgmn) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q99538

Item Number: CM0006974



Introduction

Legumain (LGMN) rabbit monoclonal antibody for Western blot, IF, IHC, ELISA. Reacts with human, mouse, rat legumain. Research use in antigen presentation and lysosomal degradation studies.

[Learn More](#)

Parameter	Specification
Product Name	Legumain (LGMN) Rabbit mAb
Remarks/Alias	LGMN; AEP; LGMN1; PRSC1; legumain; Legumain (LGMN)
Species	Human
Gene ID	5641
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-288 of human Legumain (LGMN). (NP_005597.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99538
Protein Weight	42kDa/49kDa
Shipping	Ice bag