


Anti-Legumain (Lgmn) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99538

Item Number: CM0013814



Introduction

Rabbit polyclonal antibody against human Legumain (LGMN) for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat LGMN; 49kDa. Legumain is a lysosomal cysteine protease with roles in antigen presentation and protein degradation.

[Learn More](#)

Parameter	Specification
Product Name	Legumain (LGMN) Rabbit pAb
Remarks/Alias	AEP; LGMN1; PRSC1; Legumain (LGMN)
Species	Human
Gene ID	5641
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-288 of human Legumain (LGMN) (NP_005597.3?)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99538
Protein Weight	49kDa
Shipping	Ice bag