


Anti-Leo1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Wvc0

Item Number: CM0029049



Introduction

Rabbit polyclonal antibody targeting human LEO1 (RDL) validated for WB, IHC-P, IF/ICC, and ELISA. Detects LEO1 across human, mouse, and rat samples. Ideal for transcription regulation and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	LEO1 Rabbit pAb
Remarks/Alias	RDL; LEO1
Species	Human
GeneID (human)	123169
GeneID	123169
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-570 of human LEO1 (NP_620147.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WVC0
Protein Weight	75kDa
Shipping	Ice bag