

Anti-Leptin Rabbit Polyclonal Antibody For Ihc-P, Elisa - P41159

Item Number: CM0030851



Introduction

Rabbit polyclonal antibody targeting human Leptin (OB, OBS), validated for IHC-P and ELISA. Ideal for energy balance, obesity, and diabetes research with high specificity for human Leptin.

[Learn More](#)

Parameter	Value
Product Name	Leptin Rabbit pAb
Remarks/Alias	OB; OBS; LEPD; Leptin
Species	Human
GeneID (human)	3952
GeneID	3952
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 22-167 of human Leptin (NP_000221.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P41159
Protein Weight	19kDa
Shipping	Ice bag