

Anti-Leptin Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P41159

Item Number: CM0018594



Introduction

Leptin Rabbit Polyclonal Antibody validated for Western Blot, Immunofluorescence, and ELISA. Cross-reacts with Human, Mouse, and Rat. An essential tool for studies on obesity, diabetes, and metabolic signaling.

[Learn More](#)

Parameter	Value
Product Name	Leptin Rabbit pAb
Remarks/Alias	OB; OBS; LEPD; Leptin
Species	Human
GeneID	3952 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Leptin (NP_000221.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41159
Protein Weight	19kDa
Shipping	Ice bag