

Anti-Leptin Receptor Polyclonal Antibody For Wb, If/Icc, Elisa - P48357

Item Number: CM0030554



Introduction

Rabbit polyclonal antibody targeting human leptin receptor (LEPR, CD295). Validated for WB, IF/ICC, and ELISA with cross-reactivity to rat. Ideal for investigating leptin signaling, energy homeostasis, and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	Leptin Receptor Rabbit pAb
Remarks/Alias	OBR; OB-R; CD295; LEP-R; LEPRD; Leptin Receptor
Species	Human
Gene ID (Human)	3953
Gene ID	3953
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1070-1165 of human LEPR (NP_002294.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P48357
Protein Weight	132kDa
Shipping	Ice bag