


Anti-Letm1 Rabbit Polyclonal Antibody For Wb And Elisa - 095202

Item Number: CM0020116



Introduction

High-quality rabbit polyclonal antibody targeting the N-terminal region of human LETM1 mitochondrial exchanger. Validated for WB and ELISA, cross-reacts with mouse and rat. Supports mitochondrial biology and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	LETM1 Rabbit pAb
Remarks/Alias	Mdm38; CONDMIM; SLC55A1; LETM1
Species	Human
Gene ID (Human)	3954
Gene ID	3954
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human LETM1 (NP_036450.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095202
Protein Weight	83kDa
Shipping	Ice bag