

Anti-Letmd1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q6P1Q0

Item Number: CM0024278

Introduction

Rabbit polyclonal antibody targeting human LETMD1, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial function studies and cancer research.

[Learn More](#)

Parameter	Value
Product Name	LETMD1 Rabbit pAb
Remarks/Alias	HCCR; HCCR1; HCCR2; HCCR-1; HCCR-2; HCRR-2; SLC55A3; 1110019O13Rik; LETMD1
Species	Human
GeneID	25875
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 161-360 of human LETMD1 (NP_056231.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P1Q0
Protein Weight	42kDa
Shipping	Ice bag