

Anti-Leukosialin (Spn/Cd43) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16150

Item Number: CM0027917

Introduction

Rabbit polyclonal antibody targeting human leukosialin (SPN/CD43), a leukocyte sialoglycoprotein. Suitable for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for T-cell function and leukocyte trafficking research.

[Learn More](#)

Parameter	Value
Product Name	SPN Rabbit pAb
Remarks/Alias	LSN; CD43; GALGP; GPL115; LEU-22; SPN
Species	Human
GeneID (Human)	6693
GeneID	6693
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 276-400 of human SPN (NP_003114.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16150
Protein Weight	40kDa
Shipping	Ice bag