

Anti-Lgals3Bp Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q08380

Item Number: CM0027638



Introduction

Rabbit polyclonal antibody targeting human LGALS3BP (90K/Mac-2-BP), validated for western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for cell adhesion, cancer, and immune defense research.

[Learn More](#)

Parameter	Value
Product Name	LGALS3BP Rabbit pAb
Remarks/Alias	90K; M2BP; gp90; CyCAP; BTBD17B; MAC-2-BP; TANGO10B; LGALS3BP
Species	Human
GeneID (human)	3959
GeneID	3959
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 306-585 of human LGALS3BP (NP_005558.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08380
Protein Weight	65kDa
Shipping	Ice bag