


Anti-Lgalsl Rabbit Polyclonal Antibody For Wb And Elisa - Q3Zcw2

Item Number: CM0024314



Introduction

Rabbit polyclonal antibody targeting human galectin-related protein (LGALS1). Validated for Western blotting and ELISA, with cross-reactivity to mouse. Suitable for studies on non-carbohydrate-binding lectins and galectin family functions.

[Learn More](#)

Parameter	Value
Product Name	LGALS1 Rabbit pAb
Remarks/Alias	GRP; HSPC159; LGALS1
Species	Human
GeneID (Human)	29094
GeneID	29094
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human LGALS1 (NP_054900.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q3ZCW2
Protein Weight	19kDa
Shipping	Ice bag