

Anti-Lgi1 Polyclonal Antibody For Wb, Elisa - 095970

Item Number: CM0026345



Introduction

Rabbit polyclonal antibody against LGI1, validated for WB and ELISA with reactivity to human, mouse, and rat. Useful for epilepsy and synaptic function studies.

[Learn More](#)

Parameter	Specification
Product Name	LGI1 Rabbit pAb
Remarks/Alias	EPT; ETL1; ADLTE; ADPAEF; ADPEAF; IB1099; EPITEMPIN; LGI1
Species	Human
Gene ID (Human)	9211
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 301-557 of human LGI1 (NP_005088.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95970
Protein Weight	64kDa
Shipping	Ice bag