

Anti-Lgr5 Polyclonal Antibody For Wb, Fc, Elisa - 075473

Item Number: CM0021326



Introduction

Rabbit polyclonal antibody targeting human LGR5 (GPR49, stem cell marker). Validated for WB, FC, ELISA, cross-reacts with mouse. Suitable for intestinal and hair follicle stem cell research. Shipped on ice.

[Learn More](#)

Parameter	Details
Product Name	LGR5 Rabbit pAb
Remarks/Alias	FEX; HG38; GPR49; GPR67; GRP49; LGR5
Species	Human
GeneID (Human)	8549
GeneID	8549
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 480-580 of human LGR5 (NP_003658.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human, Mouse
SWISS	075473
Protein Weight	99998 Da
Shipping	Ice bag