

Anti-Lgr5/Gpr49 Monoclonal Antibody For Wb And Elisa - O75473

Item Number: CM0002556

Introduction

Rabbit monoclonal antibody targeting human LGR5/GPR49, validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for stem cell and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	LGR5/GPR49 Rabbit mAb
Remarks/Alias	FEX; HG38; GPR49; GPR67; GRP49; LGR5/GPR49
Species	Human
Gene ID (Human)	8549
Gene ID	8549
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 800-907 of human LGR5/GPR49 (O75473)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O75473
Protein Weight	100kDa
Shipping	Ice bag