


Anti-Lgr5/Gpr49 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O75473

Item Number: CM0016023



Introduction

Rabbit polyclonal antibody targeting human LGR5/GPR49, suitable for WB, IHC-P, and ELISA. Shows cross-reactivity with mouse and rat. For stem cell and Wnt pathway studies.

[Learn More](#)

Parameter	Value
Product Name	LGR5/GPR49 Rabbit pAb
Remarks/Alias	FEX; HG38; GPR49; GPR67; GRP49; LGR5/GPR49
Species	Human
GeneID (Human)	8549
GeneID	8549
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 824-907 of human LGR5/GPR49 (NP_003658.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O75473
Protein Weight	100kDa
Shipping	Ice bag