

Anti-Lgr6 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Hbx8

Item Number: CM0009575

Introduction

Anti-LGR6 rabbit monoclonal antibody validated for Western blot, IF/ICC, and ELISA. Detects human and rat LGR6, a receptor for R-spondins that potentiates Wnt signaling and marks epidermal stem cells.

[Learn More](#)

Parameter	Value
Product Name	LGR6 Rabbit mAb
Remarks/Alias	GPCR; VTS20631; LGR6
Species	Human
GeneID (Human)	59352
GeneID	59352
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 801-900 of human LGR6 (Q9HBX8).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9HBX8
Protein Weight	104kDa
Shipping	Ice bag