

Anti-Lgr6 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hbx8

Item Number: CM0029728

Introduction

High-quality rabbit polyclonal antibody against human LGR6, validated for Western blot and ELISA. Reacts with human, mouse, and rat LGR6. Suitable for stem cell and Wnt signaling research.

[Learn More](#)

Parameter	Value
Product Name	LGR6 Rabbit pAb
Remarks/Alias	GPCR; VTS20631; LGR6
Species	Human
GeneID	59352
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-510 of human LGR6 (NP_001017403.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HBX8
Protein Weight	104kDa
Shipping	Ice bag