

Anti-Lhb Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P01229

Item Number: CM0009430

Introduction

Rabbit monoclonal antibody against human luteinizing hormone beta subunit (LHB). Suitable for immunofluorescence, immunohistochemistry, and ELISA. Ideal for reproductive endocrinology studies.

[Learn More](#)

Parameter	Value
Product Name	LHB Rabbit mAb
Remarks/Alias	CGB4; HH23; LSH-B; LSH-beta
Species	Human
GeneID (Human)	3972
GeneID	3972
Immunogen	Recombinant Protein corresponding to amino acids 21-141 of human LHB (NP_000885.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01229
Protein Weight	15 kDa
Shipping	Ice bag