

Anti-Lh(Cga/Lhb) Monoclonal Antibody For Wb, Ihc-P, Elisa - P01215 / P01229

Item Number: CM0006671



Introduction

CamelBio offers a rabbit monoclonal antibody targeting human LH (CGA/LHB), suitable for WB, IHC-P, and ELISA applications. Designed for IVD researchers developing reproductive hormone assays.

[Learn More](#)

Parameter	Value
Product Name	LH(CG/LHB) Rabbit mAb
Remarks/Alias	HCG; LHA; FSHA; GPA1; GPHA; TSHA; GPHA1; CG-ALPHA; LH(CG/LHB)
Species	Human
GeneID (Human)	1081/3972
GeneID	1081/3972
Immunogen	Recombinant proteins encompassing amino acids 25-116 of human CGA and LHB
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01215/P01229
Protein Weight	13KDa/15KDa
Shipping	Ice bag