

Anti-Lhcgr Polyclonal Antibody For Wb, If-P, Elisa - P22888

Item Number: CM0020050



Introduction

LHCGR Rabbit pAb for WB, IF-P, ELISA with cross-reactivity to mouse. Ideal for studies of lutropin-choriogonadotropic hormone receptor in reproductive biology and thyroid function. SWISS: P22888.

[Learn More](#)

Parameter	Value
Product Name	LHCGR Rabbit pAb
Remarks/Alias	LHCGR; HHG; LCGR; LGR2; LH/CG-R; LH/CGR; LHR; LHRHR; LSH-R; ULG5; lutropin-choriogonadotropic hormone receptor
Species	Human
Gene ID (Human)	3973
Gene ID	3973
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse
SWISS	P22888
Protein Weight	71kDa/78kDa
Shipping	Ice bag