

Anti-Lhcgr Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P22888

Item Number: CM0026186



Introduction

LHCGR Rabbit Polyclonal Antibody validated for Western blot, immunofluorescence, and ELISA. Reacts with human, mouse, rat LHCGR. Ideal for studying gonadotropin signaling and reproductive endocrinology. Recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	LHCGR Rabbit pAb
Remarks/Alias	HHG; LHR; LCGR; LGR2; ULG5; LHRHR; LSH-R; LH/CGR; LH/CG-R; LHCGR
Species	Human
GeneID (Human)	3973
GeneID	3973
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-360 of human LHCGR (NP_000224.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22888
Protein Weight	79kDa
Shipping	Ice bag