

Anti-Lhpp Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H008

Item Number: CM0014557



Introduction

Rabbit polyclonal antibody targeting human LHPP, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. LHPP is a phosphohistidine/phospholysine phosphatase involved in metabolism. Suitable for cancer and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	LHPP Rabbit pAb
Remarks/Alias	HDHD2B; LHPP
Species	Human
GeneID	64077
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human LHPP (NP_071409.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H008
Protein Weight	29kDa
Shipping	Ice bag