

Anti-Lhx1 Polyclonal Antibody For Wb And Elisa - P48742

Item Number: CM0022374

Introduction

Anti-LHX1 rabbit polyclonal antibody suitable for Western blot and ELISA applications. Recognizes human LHX1 with cross-reactivity to mouse and rat. Useful for studying transcription regulation and mesoderm development. UniProt P48742.

[Learn More](#)

Parameter	Value
Product Name	LHX1 Rabbit pAb
Remarks/Alias	LIM1; LIM-1; LHX1
Species	Human
GeneID (Human)	3975
GeneID	3975
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human LHX1 (NP_005559.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P48742
Protein Weight	45kDa
Shipping	Ice bag