

Anti-Lhx3 Monoclonal Antibody For Wb, Elisa - Q9Ubr4

Item Number: CM0003493

Introduction

Rabbit anti-LHX3 monoclonal antibody validated for WB and ELISA. Detects human, mouse, rat LHX3, a key transcription factor for pituitary and neuronal development. MW 43kDa.

[Learn More](#)

Parameter	Value
Product Name	LHX3 Rabbit mAb
Remarks/Alias	LIM3; CPHD3; M2-LHX3; LHX3
Species	Human
GeneID (Human)	8022
GeneID	8022
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human LHX3 (Q9UBR4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBR4
Protein Weight	43kDa
Shipping	Ice bag