

Anti-Lhx4 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q969G2

Item Number: CM0029512



Introduction

LHX4 Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, and ELISA in human, mouse, rat samples. It targets the LIM/homeobox protein Lhx4, a transcription factor involved in lung development.

[Learn More](#)

Parameter	Value
Product Name	LHX4 Rabbit pAb
Remarks/Alias	CPHD4; LHX4
Species	Human
Gene ID (Human)	89884
Gene ID	89884
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 121-390 of human LHX4 (NP_203129.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969G2
Protein Weight	43kDa
Shipping	Ice bag