

Anti-Lhx5 Rabbit Polyclonal Antibody For Wb And Elisa - Q9H2C1

Item Number: CM0019736

Introduction

Anti-LHX5 rabbit polyclonal antibody for WB and ELISA, reacting with human, mouse, rat. Targets LIM/homeobox protein Lhx5 (44 kDa, Q9H2C1), a neuronal differentiation regulator. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	LHX5 Rabbit pAb
Remarks/Alias	LHX5
Species	Human
GeneID (Human)	64211
GeneID	64211
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-180 of human LHX5 (NP_071758.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H2C1
Protein Weight	44kDa
Shipping	Ice bag