

Anti-Lhx6 Polyclonal Antibody For Wb, Elisa - Q9Upm6

Item Number: CM0029008



Introduction

Rabbit anti-LHX6 polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat LHX6. Ideal for studying GABAergic interneuron development and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	LHX6 Rabbit pAb
Remarks/Alias	hLHX6, LHX6.1, LHX6
Species	Human
Gene ID (Human)	26468
Gene ID	26468
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 143-392 of human LHX6 (NP_055183.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UPM6
Protein Weight	40kDa
Shipping	Ice bag