

Anti-Lhx6 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Upm6

Item Number: CM0021173

Introduction

LHX6 rabbit polyclonal antibody validated for WB and ELISA on human samples. LHX6 is a LIM-homeobox transcription factor critical for interneuron migration and cortical development. Suitable for neuroscience and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	LHX6 Rabbit pAb
Remarks/Alias	LHX6; LHX6.1; LIM homeobox 6
Species	Human
GeneID	26468
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human LHX6 (NP_001229263.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UPM6
Protein Weight	19kDa/38kDa/40kDa/41kDa/43kDa
Shipping	Ice bag