

Anti-Lhx8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q68G74

Item Number: CM0018248

Introduction

Anti-LHX8 rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Detects human LHX8, a LIM-homeobox transcription factor involved in neuronal and mesenchymal differentiation. Suitable for neuroscience and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	LHX8 Rabbit pAb
Remarks/Alias	LHX7; LHX8
Species	Human
Gene ID	431707
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human LHX8 (NP_001001933.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q68G74
Protein Weight	39kDa
Shipping	Ice bag