


Anti-Lhx9 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nq69

Item Number: CM0014848



Introduction

Rabbit polyclonal antibody targeting human LHX9, suitable for WB and ELISA. Cross-reacts with human, mouse, and rat. LHX9 is a LIM/homeobox protein involved in gonadal development.

[Learn More](#)

Parameter	Value
Product Name	LHX9 Rabbit pAb
Remarks/Alias	LHX9
Species	Human
Gene ID (Human)	56956
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human LHX9 (NP_064589.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQ69
Protein Weight	44kDa
Shipping	Ice bag