

Anti-Lif Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P15018

Item Number: CM0016590



Introduction

High-quality rabbit polyclonal antibody against human LIF, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat LIF (22kDa). Ideal for stem cell biology and cytokine signaling research.

[Learn More](#)

Parameter	Value
Product Name	LIF Rabbit pAb
Remarks/Alias	CDF; DIA; HILDA; MLPLI; LIF
Species	Human
GeneID (Human)	3976
GeneID	3976
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15018
Protein Weight	22kDa
Shipping	Ice bag