


Anti-Lifr Rabbit Polyclonal Antibody For Wb, Elisa - P42702

Item Number: CM0017879



Introduction

Rabbit polyclonal antibody against human LIFR, suitable for Western blotting and ELISA. Cross-reacts with mouse. Targets the leukemia inhibitory factor receptor, a signal-transducing subunit involved in cytokine receptor complexes.

[Learn More](#)

Parameter	Value
Product Name	LIFR Rabbit pAb
Remarks/Alias	SWS; SJS2; STWS; CD118; LIF-R; LIFR
Species	Human
Gene ID (Human)	3977
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human LIFR (NP_002301.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P42702
Protein Weight	124kDa
Shipping	Ice bag